

Surface Tension Measurements by Detachment and Partial Immersion Methods

Teng-Jui Lin

Department of Chemical Engineering, University of Washington

Surface and Colloid Science

Surface tension can be measured from force methods

- Detachment method (du Noüy ring)
 - Force required to pull a solid completely through a fluid interface
- Partial immersion method (Wilhelmy slide)
 - Force required to maintain the position of a solid which penetrates a fluid interface
- Force balance
 - Down = Weight - Buoyancy + Surf. Tension


$$F_{\downarrow} = F_g - F_b + F_{\sigma}$$

$$F_{\downarrow} = mg - \rho g V_{\text{disp}} + P\sigma$$


Assumes: uniform σ , fully wetted $\theta = 0$

du Noüy ring



Overall view



Side view



Wilhelmy slide



du Noüy ring can be used to measure σ by the detachment method

- Force balance

$$\begin{aligned}F_{\downarrow} &= mg - \rho g V_{\text{disp}} + P\sigma \\&= mg - 0 + 2\pi[(R + r) + (R - r)]\sigma \\&= mg + 4\pi R\sigma\end{aligned}$$

- Uncorrected surface tension σ^*



$$\sigma^* = \frac{F_{\downarrow} - mg}{4\pi R}$$

- Correction factor F [cgs unit]

$$F_{\downarrow} = mg + \frac{4\pi R\sigma}{F}$$

$$F = 0.725 + \sqrt{\frac{0.01425\sigma^*}{(2\pi R)^2(\rho_l - \rho_{\text{air}})}} + 0.04534 - \frac{1.679}{R/r}$$

du Noüy ring






Wilhelmy slide can be used to measure σ by the partial immersion method

- Force balance ($d = 0, t \ll w$)

$$\begin{aligned} F_{\downarrow} &= mg - \rho g V_{\text{disp}} + P\sigma \\ &= mg - \rho g t w d + 2(w + t)\sigma \\ &= mg + 2w\sigma \end{aligned}$$

Wilhelmy slide



Tensiometer is used to measure force on du Noüy rings and Wilhelmy slides

