

Surface Tension Measurement and Adsorption Isotherm Determination Using the Drop Weight Method

Teng-Jui Lin

Department of Chemical Engineering, University of Washington

Surface and Colloid Science

Surface tension calculation from the drop weight method

- Tate (1864)

$$W = 2\pi r\sigma$$

- Harkins and Brown (1919)

$$W = 2\pi r\sigma f = \frac{r\sigma}{F}$$

$$F \equiv \frac{1}{2}\pi f = \text{function of } (r, V)$$

- Heertjes et al. (1971)

$$F = 0.14782 + 0.27896 \left(\frac{r}{V^{1/3}} \right) - 0.1662 \left(\frac{r}{V^{1/3}} \right)^2$$

$$\text{Constraint: } \left(\frac{r}{V^{1/3}} \right) \in (0.3, 1.2)$$

- Surface tension

$$\sigma = \frac{V|\rho_2 - \rho_1|gF}{r}$$

$$\boxed{\sigma = \frac{mgF}{r}}$$

$$= \frac{\rho_2 \equiv \rho}{V} g F = \frac{m = \rho V}{r} g F$$

$$* F(r, V^{1/3})$$

need ρ

measure:
droplet's
know: ρ, V, r

replaced by m





$$W = F_{\text{net}} = |\rho_2 - \rho_1| V g$$

$$F_{\text{buoy}} = \rho_2 V g$$

$$F_g = m g = \rho_1 V g$$

$$F_g = W$$

$$\sigma = \frac{W F}{r} = \frac{|\rho_2 - \rho_1| V g}{r}$$




Szyszkowski equation describes surface tension of binary aqueous solutions

- Szyszkowski equation : $\sigma(m_2)$ ← molality (mol/kg) of the solute

$$\sigma = \sigma_0 - RTB \ln \left(1 + \frac{m_2}{a} \right)$$

↑ ↑ ↑ ↑
Surface " ideal absolute
tension of gas temp.
 pure const
 solvent
 $R = 8.314 \text{ J/mol K}$

a, B , empirical const.



Adsorption isotherm is modeled by Gibbs adsorption equation

- Gibbs adsorption equation: ideal solution

relative adsorption

$$\Gamma_{2,1} = - \frac{x_2}{RT} \frac{d\sigma}{dx_2}$$

— surface tension
— mol frac of solute

solute — solvent

- Gibbs adsorption equation: ideal dilute solution

$$\Gamma_{2,1} = - \frac{m_2}{RT} \frac{d\sigma}{dm_2}$$


— molality (mol/kg)

- General Gibbs adsorption equation

$$\Gamma_{2,1} = - \frac{m_2}{RT} \frac{d\sigma}{dm_2} \left[1 + \frac{d \ln \gamma_2^H}{d \ln m_2} \right]$$

— activity coeff based on
if $\gamma_2^H \neq f(m_2)$ Henry's law
 $\frac{d \ln \gamma_2^H}{d \ln m_2} \rightarrow 0$


alcohol aq soln



$\Gamma_{2,1} > 0 \Rightarrow$ interfacial layer is enriched of comp 2 than the bulk phase

$$\Gamma_{2,1}(m_2)$$

Objective: Determining concentration dependence of surface tension



- Surface tension of water and aqueous solutions

- $\sigma = \frac{mg}{r} F(r, V)$

- $F = 0.14782 + 0.27896 \left(\frac{r}{V^{1/3}} \right) - 0.1662 \left(\frac{r}{V^{1/3}} \right)^2$

- Surface tension of aqueous n-butanol solution vs. concentration


- Determine Szyszkowski parameters a and B

- $\sigma = \sigma_0 - RTB \ln \left(1 + \frac{m_2}{a} \right)$

$$c = \frac{m}{V}$$

$$V = \frac{m}{\rho}$$

measure σ, σ_0



known R, m_2 measure T, σ_0

calculate $\sigma(m_2)$

Objective: Constructing adsorption isotherms

$\Gamma_{2,1}$


- Relative adsorption of n-butanol at air-water interface vs. m_2 concentration (adsorption isotherm)
- Ideal dilute Gibbs adsorption equation

$$\checkmark \Gamma_{2,1} = -\frac{m_2}{RT} \frac{d\sigma}{dm_2}$$

- General Gibbs adsorption equation

$$\checkmark \Gamma_{2,1} = -\frac{m_2}{RT} \left[1 + \frac{d \ln \gamma_2^H}{d \ln m_2} \right] \frac{d\sigma}{dm_2}$$

isothermal condition
 $\Delta T = 0$



Molality m_2	$\ln(m_2)$	Activity Coefficient γ_2^H	$\ln(\gamma_2^H)$
0.003		0.9971	
0.006	A	0.9942	B
0.010		0.9906	
0.020	C	0.9823	D
0.030		0.9753	
0.040		0.9691	D-B
0.050		0.9638	C-A
0.070		0.9546	at
0.100		0.9433	$m_2 = 0.01$
0.150		0.9276	
0.200		0.9161	
0.250		0.9058	


Measure Known Calculating

T m_2, R $\Gamma_{2,1}$

$$\frac{d\sigma}{dm_2} = \frac{m_2 | \sigma}{\begin{array}{c|c} 1 & 1 \\ 2 & 2 \\ 3 & 3 \end{array}} \leftarrow \frac{\sigma_2 - \sigma_1}{m_{2,2} - m_{2,1}}$$

Objective: Determining factors affecting surface tension

- Surface tension of water vs. nozzle size
- Surface tension of water vs. drop formation time
- Surface tension of 0.01 mM Triton X-100 surfactant solution vs. drop formation time
- Compare surface tension measured by various techniques
 - Drop weight, du Noüy ring, Wilhelmy slide, sessile/pendant drop



- "actual" surface tension - "truth" thermo property - doesn't dep on measurement method
- \hookrightarrow apparent surface tension - what we measured - DOES dep on measurement method
 r, t