

# **Determination of Zeta Potential by Microelectrophoresis**


Teng-Jui Lin

Department of Chemical Engineering, University of Washington

**Surface and Colloid Science**



# Zeta potential is the electrical potential at the slip plane of the electric double layer

- $\psi$  - electrical potential difference between dispersing medium and ...
  - $\psi_\infty \equiv 0$  - the dispersing medium
  - $\psi_0$  - the surface of colloidal particle
    - True surface potential
  - $\psi_\delta$  - the outer first layer of counterions
    - Effective (Stern) surface potential
  - $\zeta$  - the slip plane (medium velocity = 0)
    - Zeta potential, electrokinetic potential
- The slip plane may be slightly further out into the solution than Stern layer
  - $\psi_\delta \approx \zeta$  ✓



# Electro-osmosis and electrophoresis give the same zeta potential

- **Electro-osmosis** - diffuse layer of ions beside a charged immobile surface move under  $E$  field, which sets the liquid into motion by the action of viscosity
  - *Liquid* moves, *solid* at rest
- **Electrophoresis** - diffuse layer of ions beside a charged particle surface move under  $E$  field, which sets the particle into motion
  - *Solid* moves, *liquid* at rest



# Electrophoretic mobility allows determination of zeta potential

- **Electrophoretic mobility** -  $u_E [(\mu\text{m/s})/(\text{V/cm})]$

$$\circ u_E = \frac{V_p}{E_x} = \begin{cases} \frac{\varepsilon\varepsilon_0\zeta}{\mu} & (\kappa a > 200) \\ \frac{2}{3} \frac{\varepsilon\varepsilon_0\zeta}{\mu} & (\kappa a < 0.1) \end{cases}$$

Helmholtz-Smoluchowski limit  
Huckel limit

- $u_E$  - Electrophoretic mobility
- $V_p$  - Particle velocity Measure
- $E_x$  - Electric field strength Known
- $\varepsilon$  - Dielectric constant of the medium Known - H<sub>2</sub>O
- $\varepsilon_0$  - Permittivity of free space Known
- $\zeta$  - Zeta potential Want
- $\mu$  - Viscosity of the medium Known - H<sub>2</sub>O
- $\kappa^{-1}$  - Debye length
- $a$  - Particle radius


measure  
↓

$$\zeta = \frac{V_p \mu}{E_x \varepsilon \varepsilon_0}$$

✓   ✓   ✓   ✓

# Point of zero charge and isoelectric point define pH at which potentials are zero

---


- **Potential determining ions** - ions whose concentration determines surface potential
  - Crystalline solid - lattice ions
  - Oxides -  $\text{H}_3\text{O}^+$ ,  $\text{OH}^-$  (pH)
- **Point of zero charge (PZC)** - pH at which  $\psi_0 = 0$ 
  - $\text{pH} < \text{PZC}$ :  $\psi_0 > 0$
  - $\text{pH} > \text{PZC}$ :  $\psi_0 < 0$
- **Isoelectric point (IEP, pl)** - pH at which  $\zeta \approx \psi_\delta = 0$ 
  - $u_E = 0 \Rightarrow V_p = 0$  

# Darkfield illumination microscopy visualizes colloidal particles under $E$ field


- $u_E = \frac{V_p}{E_x} = \frac{\epsilon\epsilon_0\zeta}{\mu} \quad (\kappa a > 200)$

- $\zeta = \frac{\mu V_p}{\epsilon\epsilon_0 E_x}$  particle velocity

- Electrophoresis has solid particles moving when liquid is at rest ( $v_x = 0$ )



- $v_x = 0$  when  $y = \begin{cases} 0.205w & \checkmark \\ 0.795w & \checkmark \end{cases}$



# Laser Doppler electrophoresis determines zeta potential with more sensitivity

