

# **BET (Brunauer-Emmett-Teller) Method for Surface Area Determination**

Teng-Jui Lin

Department of Chemical Engineering, University of Washington

**Surface and Colloid Science**

# BET isotherms empirically fits adsorption isotherm data

- BET method - procedure for determining surface area of finely divided solid by measurements of low-temperature gas adsorption

- BET isotherm

$$\frac{V}{V_m} = \frac{C(p/p_0)}{(1 - p/p_0)[1 + (C - 1)(p/p_0)]}$$

*vol of gas adsorbed* (for  $V$ )  
*vol of gas if monolayer* (for  $V_m$ )  
*const pressure of adsorbate in gas* (for  $p/p_0$ )  
*sat press.* (for  $p_0$ )

- Linear form of BET isotherm

$$\left[ \frac{(p/p_0)}{V(1 - p/p_0)} \right] = \frac{1}{CV_m} + \frac{(C - 1)(p/p_0)}{CV_m}$$



*y* (for the first term)  
*x* (for the second term)

- Single point method:  $C \gg 1$ , so intercept  $\approx 0$

Volume:  $V_m = V(1 - p/p_0)$

Mass:  $\underline{X}_m = X(1 - p/p_0)$

$$\leftarrow \frac{1}{V(1 - p/p_0)} = \frac{1}{V_m}$$



*monolayer (Langmuir)*

# BET method allows surface area determination

- Volume of adsorbate desorbed by sample

$$V_{\text{des}} = V_c \frac{A}{A_c}$$

*calibration* (pointing to  $A$ )  
*sample integrator* (pointing to  $A$ )  
*calibration integrator* (pointing to  $A_c$ )  
*measurable* (bracketed around the fraction)

- Mass of adsorbate desorbed *ideal gas law*

$$X = \frac{p_a V_{\text{des}}}{RT} M_a = \frac{p_a V_c (A/A_c)}{RT} M_a$$

- Mass of adsorbate desorbed when solid is covered by a single adsorbate monolayer

$$X_m = \frac{p_a V_c (A/A_c)}{RT} M_a (1 - p/p_0)$$

$$X_m = \sum (1 - p/p_0)$$

*←* (arrow from the previous equation)

- Total surface area

$$S_t = X_m \frac{N_0}{M_a} A_{cs}$$

*measured* (circled around  $S_t$ )  
*Arg #* (pointing to  $N_0$ )  
*# molecules* (pointing to  $N_0$ )  
*Cross-sec area of molecule* (pointing to  $A_{cs}$ )

- Specific surface area

$$\Sigma = \frac{S_t}{m}$$

*← mass of sample* (pointing to  $m$ )

# Initial configuration of FlowSorb Surface Area Analyzer

- Initial configuration
  - Power: OFF
  - Toggle valve: open
  - Sample holders are installed at DEGAS, TEST, and COLD TRAP locations
  - Insulating container is installed at sample holder for COLD TRAP



# Startup procedure of FlowSorb Surface Area Analyzer

---

- Startup
  - Open inert gas flow (30% nitrogen, 70% helium), wait 5 min
  - Adjust flow meter to the calibration mark
  - Pour liquid nitrogen into an insulating container, and pour the liquid nitrogen to the insulating container at COLD TRAP, wait for 10 min
  - Power switch: ON, wait for 30 min

# Preparing activated carbon sample for BET measurement

---

- Sample preparation
  - Weigh 0.02 g activated carbon in the sample holder
  - Place the sample holder on DEGAS and heat with a heating mantle, wait for 15 min (concurrent with start up)

# Calibrating FlowSorb Surface Area Analyzer

---

- Calibration
  - Flush gas syringe with evaporated nitrogen gas above the liquid nitrogen at DEGAS
  - Fill the gas syringe with 1 mL nitrogen gas
  - Wipe needle tip free of frost, set aside needle to equilibrate to room temperature
  - Zero the instrument display with course and fine zeros. Switch to surface area (SA), and clear SA display
  - Insert needle at INJECT, inject nitrogen gas at moderate rate, and withdraw
  - Wait until reading is stabilized, then calibrate the instrument to **SA = 2.84**
    - Confirm calibration by repeated injection, if necessary

# Measuring surface area of activated carbon with the BET method (adsorption)

---

- Adsorption measurement
  - Exchange sample holders between DEGAS and TEST, so the sample is at TEST
  - Wait until readings stabilize and clear SA display
  - Pour liquid nitrogen into an insulating container
  - Place the insulating container at TEST by immersing the sample holder with liquid nitrogen. Secure the container by flipping on the container holder
    - Nitrogen gas starts to adsorb to activated carbon
  - Wait until reading is stabilized, **record the value of SA** = adsorption surface area



# Measuring surface area of activated carbon with the BET method (desorption)

---

- Desorption measurement
  - Remove the insulating container by putting down the container holder
    - Nitrogen gas starts to desorb from activated carbon
  - Wait until reading is stabilized, **record the value of SA** = desorption surface area

# Shutdown procedure of FlowSorb Surface Area Analyzer

---

- Shutdown
  - Power: OFF
  - Remove the activated carbon from the sample holder, and place the sample holder back on TEST
  - Turn off the inert gas flow

# Calculating the specific surface area of activated carbon

---

- Measured

$S_t$ : Adsorption SA, desorption SA

$m$ : Mass of activated carbon

- Calculated

- $\Sigma = \frac{S_t}{m}$  - specific SA for adsorption and desorption

Usually use desorption data to report  $\Sigma$