

# **Interfacial Tension Measurement and Adsorption Isotherm Determination Using the Inverted Drop Weight Method**

Teng-Jui Lin

Department of Chemical Engineering, University of Washington

**Surface and Colloid Science**

# Surface tension calculation from the drop weight method

- Tate (1864)

$$W = 2\pi r\sigma \quad \leftarrow$$

- Harkins and Brown (1919)

$$W = 2\pi r\sigma f = \frac{r\sigma}{F}$$

$$F \equiv \frac{1}{2}\pi f$$

$$f = \frac{\text{W of liq breaks away}}{\text{W of pendant drop}}$$

$$\sigma = \frac{WF}{r} = \frac{|\rho_1 - \rho_2| g V \cdot F}{r}$$

$$\sigma P = \sigma 2\pi r$$



$$= \rho_1 g V$$

$$F_{\text{buoy}}$$

$$F_g = W$$

$$W = mg$$

$$= \rho_2 V g$$



- Heertjes et al. (1971)

$$F = 0.14782 + 0.27896 \left( \frac{r}{V^{1/3}} \right) - 0.1662 \left( \frac{r}{V^{1/3}} \right)^2 = \underline{f(r, V)}$$



Constraint:  $\left( \frac{r}{V^{1/3}} \right) \in (0.3, 1.2)$

- Surface tension


$$\sigma = \frac{V|\rho_2 - \rho_1|gF}{r}$$

if  $\rho_1 \rightarrow 0$ ,  $\rho_2 \equiv \rho$

$$= \frac{V\rho g F}{r} = \frac{mg F}{r}$$



# Interfacial tension between liquids is measured with inverted drop weight method



- Interfacial tension  $F = f(r, V)$ 
  - $\sigma = \frac{V|\rho_2 - \rho_1|gF}{r}$
  - $F = 0.14782 + 0.27896 \left( \frac{r}{V^{1/3}} \right) - 0.1662 \left( \frac{r}{V^{1/3}} \right)^2$





# Szyszkowski equation describes surface tension of binary aqueous solutions

- Szyszkowski equation

$$\sigma = \sigma_0 - RTB \ln \left( 1 + \frac{C_2}{a} \right)$$

$\sigma$ : surface tension  
 $\sigma_0$ : " of pure solvent  
 $R$ : ideal gas const  
 $T$ : abs temp  
 $B$ : empirical const  
 $C_2$ : concentration (molarity) of solute  
 $a$ : empirical const

$R = 8.314 \text{ J/mol K}$



# Adsorption isotherm is modeled by Gibbs adsorption equation

- Gibbs adsorption equation: ideal dilute solution

relative adsorption  
of solute  
with respect to solvent

$$\Gamma_{2,1} = -\frac{C_2}{RT} \left( \frac{d\sigma}{dC_2} \right)$$

- Finite difference method

$$\frac{d\sigma}{dC_2} = \frac{\Delta\sigma}{\Delta C_2} = \frac{\sigma_2 - \sigma_1}{C_{2,2} - C_{2,1}}$$

- Szyszkowski equation

$$\sigma = \sigma_0 - RTB \ln \left( 1 + \frac{C_2}{a} \right)$$

$$\frac{d\sigma}{dC_2} = \frac{d}{dC_2} \left( \sigma_0 - RTB \ln \left( 1 + \frac{C_2}{a} \right) \right)$$

$$= -RTB \frac{1/a}{1 + \frac{C_2}{a}}$$

$$\hookrightarrow \Gamma_{2,1} = -\frac{C_2}{RT} \frac{d\sigma}{dC_2} = +\frac{C_2}{RT} \left( +RTB \frac{1}{a} \frac{1}{1 + \frac{C_2}{a}} \right) = \frac{C_2 B}{a(1 + \frac{C_2}{a})} = \frac{C_2 B}{a + C_2}$$

data - noisy

$C_2$	$\sigma$	$\frac{d\sigma}{dC_2}$
1	1	
2	2	
3	3	
4	4	

$\frac{3-1}{3-1} = 1$

