

Centrifugation

Teng-Jui Lin

Department of Chemical Engineering, University of Washington

Separation Processes

Sedimentation is driven by gravity and follows Stoke's law at low velocity

- Sedimentation separates solid particles dispersed in liquid
- **Stoke's law** describes flow around a sphere at low velocity ($Re < 1$)

- $$v_{\infty} = \frac{(\rho_p - \rho_f)}{18\mu} D_p^2 g$$
 for $Re = \frac{D_p v_{\infty} \rho_f}{\mu} < 1$

- v_{∞} - terminal velocity
 - ρ_p - particle density
 - ρ_f - fluid density
 - μ - fluid viscosity
 - D_p - particle diameter
- Sedimentation is slow... How can we speed it up?

Centrifugation is driven by centrifugal force and also follows Stoke's law

- Centrifugation separates solid particles dispersed in liquid *faster*
- **Stoke's law** describes flow around a sphere at low velocity ($Re < 1$)

$$\circ \quad v_{\infty} = \frac{(\rho_p - \rho_f)}{18\mu} D_p^2 \omega^2 r \quad \text{for} \quad Re = \frac{D_p v_{\infty} \rho_f}{\mu} < 1$$

- ω - angular velocity of centrifuge bowl
- r - centrifuge bowl radius



- **G-force** normalizes centrifugal driving force by gravitational driving force

$$\circ \quad Z = \frac{\omega^2 r}{g}$$

Sigma factor compares performance of continuous centrifuges of the same type

- **Sigma factor** - effective area of a continuous centrifuge

- $\Sigma = \frac{\dot{V}}{2v_{\infty}}$

- **Ex.** At small scale, cells can be centrifuged at $\Sigma_1 = 200$ and $\dot{V}_1 = 15$ mL/min. At large scale, cells can be centrifuged at $\Sigma_2 = 9000$ and $\dot{V}_2 = 700$ mL/min at the same speed. The densities and viscosity are unchanged. Quantify the changed physical property.